














DIALAB Produktion und Vertrieb von chemisch-technischen Produkten
und Laborinstrumenten Gesellschaft m.b.H.
Hondastrasse, Objekt M55, IZ NOE-Sued, 2351 Wiener Neudorf, Austria,
Phone: +43 (0) 2236 660910-0, Fax: +43 (0) 2236 660910-30, e-mail:
office@dialab.at

Instruction for Use

RINSE®


Explanation of symbols

	Reference Catalog
	Expiration date
	Lot number
	Manufacturer
	Consult the instructions for use
	Compliance with European Directive 98/79 / EC
	When followed by "WARNING", this symbol means "Attention! Failure to do so may lead to erroneous examination results and / or degradation of the device." "
	In vitro diagnostic medical device
	Protect from the sun
	Keep away from moisture
	Temperature limits (lower and upper)

Introduction

This document provides the information necessary for the use of RINSE® with automatic MS4-type hematology analyzers manufactured by MELET SCHLOESING Laboratoires.

This document is provided in electronic format (e-ifu), using the following URL:
<https://www.dialab.at/en/news/info/>

 **WARNING** This user manual contains important information regarding the safety of using RINSE®. Please read this document as a whole before using the product.

Identification of the legal manufacturer



DIALAB Produktion und Vertrieb von chemisch-technischen Produkten und Laborinstrumenten
Gesellschaft m.b.H.
Hondastrasse, Objekt M55, IZ NOE-Sued
2351 Wiener Neudorf, Austria
Tél : +43 (0) 2236 660910-0
Fax : +43 (0) 2236 660910-30

Product designation	CE IVD	product reference
RINSE®		RMSRINC2L0 RMSRINC500

Description
RINSE® is a rinsing agent suitable only for MS4 automatic blood cell analyzers of the brand MELET SCHLOESING Laboratories. This reagent allows the rinsing of the apparatus and the desalination of the parts in contact with the diluent.
Kit Contents
RINSE® is combined in different capacities: a) RMSRINC2L0 bottle of 2 l b) RMSRINC500 500 ml bottle.
Indication for Use
RINSE® is an in vitro diagnostic device intended for use by healthcare professionals and only in combination with the conductimetric hematology devices of the MS4 range manufactured by MELET SCHLOESING Laboratories.
Intended Use
RINSE® reagent is a rinsing agent designed to keep the device clean.
Composition
RINSE® is composed of the following elements: <ul style="list-style-type: none"> • • Conservative (<1%)
Storage conditions - Shelf life - Stability
RINSE® must never be subjected to temperatures below 15 ° C and above 25 ° C. RINSE® can be stored for 120 days (16 weeks) after opening. RINSE® can be stored for 3 years (see expiration date on the product label) since the date of manufacture.
Installation
RINSE® is ready to use and requires no preparation. A detailed description of the reagent installation / replacement is available in the relevant analyzer user guide.

Collection and preparation of the sample

When collecting the blood sample, avoid any violent aspiration to avoid the phenomenon of hemolysis that could influence the results of the hematological analysis. Also be sure to limit the sampling time to avoid micro-coagulation problems.

The blood sample to be analyzed should be collected in a collection tube containing EDTA K3 anticoagulant. The use of the sampling tube must be carried out according to the instructions of the supplier.

Gentle and prolonged homogenization of the blood / anticoagulant mixture is essential before any analysis, refer to the supplier's instructions.

It is recommended to carry out the analysis no earlier than 30 minutes and no later than 8 hours after sampling.

Test principle

The combination of the RINSE ® reagent with the MS4 meter-analyzers from MELET SCHLOESING LABORATOIRES makes it possible to rinse the parts in contact with a salt reactant such as the diluent.

Disposal of the product

Respect the protocol applied by your laboratory for the neutralization and elimination of waste

Warnings and precautions

Precautions:

- Use only for in vitro diagnostics.
- Please read the Safety Data Sheet (MSDS) before use.
- Please follow all pre-analytical steps in the laboratory.
- Handle reagents with care to prevent the formation of bubbles.
- Do not use directly after transport or directly after handling. Ensure that the operating temperature conditions are reached.

Warnings:

- Reagents can cause irritation of the eyes, skin and mucous membranes.
- In case of contact, rinse with plenty of water and seek immediate medical attention. In case of accidental ingestion, call a doctor immediately!
- Avoid particles or bacteria coming into contact with the reagent.
- Do not use reagents beyond the expiry date or beyond the open bottle time.
- Place the reagents at the main unit of the device.
- Do not transfer / mix reagents whether they are of the same type or batch.
- Do not reuse an empty container as this may distort the results or damage the machine.
- Do not use the product in case of deterioration of the protective packaging.
- If a reagent has been accidentally spilled, wipe it off immediately with a damp cloth.
- Do not use the product if there is any sign of deterioration (turbidity, color change, etc.)